

## Dexnyl® PEEK PF 8800FP

### Technical test report

| Test Items            | Test Method | Unit      | Act.Value |
|-----------------------|-------------|-----------|-----------|
| Tm                    | ISO11357    | °C        | 344       |
| Tg                    | ISO11357    | °C        | 149       |
| D50) Particle Size    | ISO13320    | D50/μm    | 26.31     |
| Bulk Density          | GB/T23771   | g/100ml   | 29        |
| Water Content         | PFL-TM-04   | %         | 0.20      |
| Melt Volume Flow Rate | ISO 1133    | cm3/10min | 81.3      |

Statement: This certificate has been prepared with care and the best of our knowledge as part of the testing system of our supplier. However, it does not relieve our customers of their obligation to inspect the goods upon receipt or establish any warranties to third parties, and it is only a reference to the general basic performance information of the product. The document issued electronically, valid without signature.

#### Distribution by:

BIEGLO GmbH  
Bahrenfelder Straße 242  
22765 Hamburg  
+49 40 4011 30000

[granules@bieglo.com](mailto:granules@bieglo.com)  
[www.bieglo.com](http://www.bieglo.com)  
[www.peek-shop.com](http://www.peek-shop.com)